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### 3.1 Preparation of the ion microsensor

Ca<sup>2+</sup>- microsensor were prepared as described previously<sup>[4]</sup>. Pre-pulled and silanized microsensor ( $\Phi 5 \pm 1 \mu\text{m}$ , XY-CGQ-01<sup>①</sup>, YoungerUSA) were first filled with a backfilling solution (100 mM CaCl<sub>2</sub><sup>②</sup>) to a length of approximately 1.0 cm from the tip.

The micropipettes were front filled with 15-50<sup>③</sup>  $\mu\text{m}$  columns of selective liquid ion-exchange cocktails (Ca<sup>2+</sup> LIX, XY-SJ-Ca<sup>④</sup>, YoungerUSA). An Ag/AgCl wire microsensor holder YG003-Y11, (YoungerUSA) was inserted in the back of the microsensor to make electrical contact with the electrolyte solution. YG003-Y11, (YoungerUSA) was used as the reference microsensor.

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### 3.2 Preparation of oxygen microsensor

The Pt/Ir polarographic oxygen microsensor (tip diameter 2-4  $\mu\text{m}$ , XY-DJ-501<sup>①</sup>, YoungerUSA) were used for detecting dissolved oxygen in -750mV<sup>①</sup> polarization voltage as described earlier<sup>[4]</sup>. YG003-Y11, (YoungerUSA) was used as the reference microsensor to complete the circuit.

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### 4.1 Calibration of the ion microsensor

Prior to the flux measurement, the microsensor were calibrated with cultural media having different concentrations of Ca<sup>2+</sup>, 0.1 mM and 0.01 mM<sup>①</sup> respectively. Only microsensor with a Nernstian slope >22<sup>②</sup> mV/decade were used in this study. The same microsensor were calibrated again according to the same procedure and standards after each test. Data was discarded if the post-test calibrations failed.

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<p><b>4.2 Calibration of oxygen microsensor</b></p> <p>The oxygen microsensor was calibrated with nitrogen- saturated and common cultural media.</p>	<p>参考文献 4</p> <p>以 O<sub>2</sub> 为例，其他的分子（如 H<sub>2</sub>O<sub>2</sub>，IAA）校正时分别采用的是含有该分子不同浓度的测试液。</p>
<p><b>5.Sampling rules</b></p> $J = -D_0 \cdot (dc/dx)$ <p><i>dx</i> is the distance the microsensor moved repeatedly from one point to another perpendicular to the surfaces of samples at a frequency of ca. 0.3Hz. Moreover, <i>dx</i> is usually between 5 to 35μm, and we used 30<sup>①</sup> μm throughout the experiments.</p>	<p>参考文献 8</p> <p>①测试时先在近样品处测得一个电位，然后远离样品 <i>dx</i> 再测得一个电位，二者作差得到实验的原始数据—电位差。<i>dx</i> 即为采样规则，一般为 5-35μm，具体的采样规则请以您的实验为准，分子的采样规则与离子相同。</p>
<p><b>6.Flux data analysis</b></p> <p>The micro-volt differences (ΔμV) were then exported as raw data before they were imported and converted into net Ca<sup>2+</sup> fluxes by using the <a href="#">JCal V3.3</a> ( a free MS Excel spreadsheet, <a href="#">youngerusa.com</a> or <a href="#">xbi.org</a>).</p>	<p>参考文献 3</p> <p>JCal V3.3 包括所有离子/分子的数据换算方法，只需要在其中找到对应的离子/分子进行计算即可。具体的流速换算表版本请以您使用的版本号为准。</p>

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1. 印莉萍等.非损伤性扫描离子选择电极技术及其在高等植物研究中的应用.自然科学进展,2006,16:262-266.
2. Li M, *et al.* Glycine betaine-mediated potentiation of HSP gene expression involves calcium signaling pathways in tobacco exposed to NaCl stress. *Physiologia Plantarum*, 2014, 150(1): 63-75.
3. Bai L, *et al.* A Receptor-Like Kinase Mediates Ammonium Homeostasis and Is Important for the Polar Growth of Root Hairs in Arabidopsis. *Plant Cell*, 2014, 26(4): 1497-1511.
4. Xu Y, *et al.* Application of non-invasive microsensing system to simultaneously measure both H<sup>+</sup> and O<sub>2</sub> fluxes around the pollen tube. *Journal of Integrative Plant Biology*, 2006,48(7): 823- 831.
5. Gilad Twig, *et al.* Real-time detection of reactive oxygen intermediates from single microglial cells. *Biological Bulletin*,2001,201:261-262.
6. McLamore ES, *et al.* Non-invasive quantification of endogenous root auxin transport using an integrated flux. *Plant Journal*.2010,63(6):1004-1016.
7. Stefano Mancuso, *et al.* Noninvasive and continuous recordings of auxin fluxes in intact root apex with a carbon nanotube-modified and self-referencing microelectrode. *Analytical biochemistry*,2005,341:344-351.
8. 孙健.胡杨响应盐胁迫与离子平衡调控信号网络研究.北京:北京林业大学,2011:26-27.
9. Chen Z, *et al.* Root Plasma Membrane Transporters Controlling K<sup>+</sup>/Na<sup>+</sup> Homeostasis in Salt-Stressed Barley. *Plant Physiology*.2007,145(4):1714-1725.

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